

# Clearview Complete HIV 1/2 Package Insert Questionnaire

Please answer the questions below.

In the right-hand column, also note the page of the package insert where you found the answer.

**We suggest that you underline or highlight the answers in your package insert for future reference.**

Question	Package Insert Page #
<b>COMPLEXITY</b>	
1. This <b>waived</b> complexity rapid HIV test can be used for the following specimen types (circle all that apply): <ul style="list-style-type: none"> <li>• Fingerstick whole blood</li> <li>• Venipuncture whole blood</li> <li>• Oral fluid</li> <li>• Plasma</li> <li>• Serum</li> </ul>	Page _____
<b>INTENDED USE</b>	
2. This test detects antibodies to (circle all that apply): <ul style="list-style-type: none"> <li>• HIV-1</li> <li>• HIV-2</li> <li>• HBV</li> <li>• HCV</li> </ul>	Page _____
<b>RESTRICTIONS</b>	
3. Test providers must use the manufacturer's instructional materials. True False	Page _____
4. Clients must receive the "Subject Information" pamphlet before specimen collection and appropriate information when test results are provided. True False	Page _____
<b>SUMMARY AND EXPLANATION OF THE TEST</b>	
5. The absence of antibody to HIV is absolute proof that an individual is free of HIV infection and incapable of transmitting the virus. True False	Page _____
<b>MATERIALS PROVIDED</b>	
6. Are kit controls and retractable lancets included in the Test Kit? True False	Page _____

<b>WARNINGS</b>	
7. This test should be performed at temperatures in the range of (circle all that apply): <ul style="list-style-type: none"> <li>• 45 – 100 degrees Fahrenheit</li> <li>• 59 – 80.6 degrees Fahrenheit</li> <li>• 18 – 30 degrees Celsius</li> <li>• 64 – 86 degrees Fahrenheit</li> </ul>	Page _____
8. Clients infected with HIV-1 and/or HIV-2 who are receiving highly active antiretroviral therapy (HAART) may produce false negative results. True False	Page _____
<b>SAFETY PRECAUTIONS</b>	
9. Which of the following statements about handling and testing <b>blood</b> specimens are correct (circle all that apply): <ul style="list-style-type: none"> <li>• Handle as if capable of transmitting infectious agents</li> <li>• Wearing disposable gloves is optional</li> <li>• Dispose of used gloves in regular trash</li> <li>• Do not eat, drink, or smoke in the area</li> </ul>	Page _____
<b>HANDLING PRECAUTIONS</b>	
10. Which of the following statements about handling instructions are correct (circle all that apply): <ul style="list-style-type: none"> <li>• Testing Devices can be re-used.</li> <li>• Do not use the test beyond the expiration date.</li> <li>• The expiration date is printed on the pouch.</li> <li>• Do not mix reagents from different kit lot numbers.</li> <li>• Use adequate lighting when reading results</li> </ul>	Page _____
<b>STORAGE INSTRUCTIONS</b>	
11. Unused Test Kits can be stored at (circle all that apply): <ul style="list-style-type: none"> <li>• 59 – 99 degrees Fahrenheit</li> <li>• 46 – 86 degrees Fahrenheit</li> <li>• 8 – 30 degrees Celsius</li> <li>• 15 – 32 degrees Celsius</li> </ul>	Page _____
<b>GENERAL TEST PREPARATION</b>	
12. When Kit contents are stored in the refrigerator, they must be brought to room temperature (64 degrees – 86 degrees Fahrenheit / 18 – 30 degrees Celsius). True False	Page _____
<b>SPECIMEN COLLECTION AND TESTING PROCEDURE</b>	
13. If there is no absorbent packet included with the Test Device: <ul style="list-style-type: none"> <li>• Report it to the manufacturer</li> <li>• Discard the Device and test with a new Pouch</li> <li>• Perform the test using the Test Device with the missing packet.</li> </ul>	Page _____
<b>QUALITY CONTROL</b>	

<p>14. A pink-purple line within the control area of the Test Device must be present in all valid tests: True False</p>	<p>Page _____</p>
<p>15. If the built-in Control Line is present, it indicates that (circle all that apply):</p> <ul style="list-style-type: none"> <li>• A specimen was added to the Vial.</li> <li>• The fluid migrated through the Test Device.</li> <li>• Test is valid whether or not the sample is reactive or non-reactive.</li> </ul>	<p>Page _____</p>
<b>TEST RESULT AND INTERPRETATION OF TEST RESULT</b>	
<p>16. A reactive test result means that (circle all that apply):</p> <ul style="list-style-type: none"> <li>• HIV-1 or HIV-2 antibodies have been detected.</li> <li>• The Control Line was absent; the Test Line was absent.</li> <li>• The Test Line is interpreted as preliminary positive.</li> </ul>	<p>Page _____</p>
<p>17. If a test is invalid (circle all that apply):</p> <ul style="list-style-type: none"> <li>• It can not be interpreted.</li> <li>• Repeat the test with a new pouch and sample.</li> <li>• Report the test as negative.</li> </ul>	<p>Page _____</p>
<b>LIMITATIONS OF THE TEST</b>	
<p>18. A non-reactive result can occur if the client has recently been exposed to HIV. True False</p>	<p>Page _____</p>
<p>19. A person who has participated in an HIV vaccine study may develop antibodies to HIV. True False</p>	<p>Page _____</p>